



VisuMax

Defining the pulse rate in refractive surgery



We make it visible.



//VisuMax
MADE BY CARL ZEISS

Its remarkable features are its precision and innovative detail.

New trends in modern corneal surgery.

With the VisuMax®, Carl Zeiss is significantly shaping the world of refractive surgery. This ground-breaking laser system employs high-performance femtosecond laser technology and is characterized by its outstanding cutting precision, unsurpassed speed and gentle treatment technique. The VisuMax is thus the ideal platform for therapeutic and refractive applications of cutting-edge corneal surgery, including Flap, Keratoplasty, Incision for ICR and ReLEx®.

With ReLEx, VisuMax heralds a paradigm shift in refractive surgery: minimally invasive laser vision correction. It completes the unparalleled range of innovative surgical possibilities and creates the ideal preconditions for tapping into new patient groups.

VisuMax is the logical enhancement of the product range for refractive surgeons – and marks another step into the future of corneal surgery, the progress of which Carl Zeiss has been shaping for more than 20 years.

VisuMax applications

Precision in all facets

The VisuMax® is extraordinary in all major applications of state-of-the-art corneal surgery. The innovative femtosecond laser system brings together perfectly coordinated components to provide maximum cutting precision, efficiency, predictability and comfort.

ReLEx

With the minimally invasive refractive procedure, ReLEx® smile, VisuMax enables a new, flapless technique. In a single step, the femtosecond laser creates the refractive lenticule, and the access through which the lenticule is extracted. Without ablation and without a flap.

Flap

In Femto-LASIK and Laser Blended Vision for treating presbyopic patients, the VisuMax stands out as a flap cutter. It provides predictable flap thickness and adjustable geometries.

Keratoplasty

VisuMax offers a broad spectrum of corneal transplant procedures, including lamellar and penetrating keratoplasty. High-precision cutting quality and rapid incision progress enable the preparation of precision corneal grafts and ideal preparation of the recipient's cornea.

Incision for ICR

The VisuMax also impresses with its femtosecond laser technology for the implantation of intracorneal rings. It even permits inclined cutting geometries and ring tunnels smaller than 360° and offers a previously unachieved degree of flexibility when defining the tunnel parameters.



VisuMax strengths

Building blocks of state-of-the-art femtosecond technology



A contact glass as ingeniously designed as the cornea

The surface of the human cornea is curved. Thus, Carl Zeiss contact glasses are curved way, too. The three different sizes available (S, M, L) ensure an optimal fit to the anatomy of the eye. The cornea is not forced into a deformed planar, non-physiological shape – and artifacts are avoided in the cutting result, as is unnecessarily high IOP.




Maximum cutting precision

High-precision optics from Carl Zeiss provide an extremely focused laser beam. The result: minimum laser pulse energy at a high pulse frequency for hitherto unsurpassed incision control – at precisely the desired depth in the cornea, even with three-dimensional, curved incisions.



Brilliant visual control

The integrated, high-quality ZEISS surgical microscope, including digital video camera for live recording of the surgical procedure ensures precise and complete control of each treatment step.



A smart unit

The ergonomically pivoting patient supporting system ensures maximum comfort. The patient's position is continuously monitored during treatment and the sturdy yet comfortable patient supporting system is automatically adjusted during surgery.


Achieving objectives with intuition

The VisuMax can be easily controlled via its touch-screen and intuitive software. An interactive wizard supports the surgeon during all steps.



Efficiency paying off

With a laser pulse frequency of 500 kHz, the VisuMax realizes short treatment times. This means more comfort for both physician and patient. In addition, the user benefits from a more efficient workflow and a higher throughput of satisfied patients.



Slit illumination for greater detail

As a universal workstation for corneal surgery, the VisuMax features integrated slit illumination for monitoring and control immediately after the respective treatment, without the patient having to switch places or change position.

ReLEx

For a new flapless treatment

With ReLEx®, the VisuMax® enables the fusion of cutting-edge femtosecond technology and precise lenticule extraction to provide minimally invasive vision correction. A refractive lenticule is created in the intact cornea and removed via a small incision. Without ablation. Without a flap. The treatment is thus a **flapless, all-femto** and **single-step** procedure.



Flapless

ReLEx is turning the world of refractive surgery on its head. In ReLEx smile, a small incision is sufficient to extract the lenticule. The minimally invasive incision also implies fewer transected nerves and significantly reduced incidence of dry-eye syndrome. The risk of infections, epithelial ingrowths or flap complications can be reduced. The smaller incisions enable better healing of the epithelium.

All-femto

With the VisuMax, ReLEx relies completely on femtosecond technology. The unique laser vision correction procedure creates the pre-calculated lenticule in the cornea with precision and predictability. Nomograms or fluence tests are not required, intraoperative ambient conditions or individual corneal characteristics have virtually no influence on the reproducibility of the lenticule. Physicians benefit from excellent predictability, particularly when correcting higher refraction values.

Single-step


With ReLEx, both the lenticule and the access incision are created in a single treatment step. Unlike Femto-LASIK, therefore, only one surgical procedure needs to be planned. There is no need for the patient to switch places. The result is more efficient workflows and shorter treatment times. For patients, the procedure is a much less stressful experience.

Outstanding results


- Small incision of less than 4 mm
- Side-cut length up to 80% shorter and
- Cap incision area up to 30% smaller than for a Femto-LASIK flap
- Implies lower incidence of dry eye syndrome and less nerve transection due to smallest of incisions with no flap
- Less risk of infections, epithelial ingrowths and flap complications
- Reproducibility of the lenticule, irrespective of individual corneal characteristics and ambient conditions
- Excellent predictability, particularly when correcting higher refraction values
- Efficient treatment process without patient having to switch places

ReLEx smile

Small Incision Lenticule Extraction



In a single step, the VisuMax creates a refractive lenticule and a small incision of less than 4 mm in the otherwise intact cornea - almost irrespective of the ambient conditions or corneal structure.



The lenticule is removed through the small incision. The disruption to the biomechanics of the cornea is minimal. No flap is cut.



The removal of the lenticule alters the shape of the cornea, thus achieving the desired refractive change.

ReLEx is not intended for sale in the United States.

Flap

For best possible conditions

For Femto-LASIK and the treatment with Laser Blended Vision, VisuMax® means above all one thing: high-precision flaps. Combined with the excimer laser MEL 80™ and the treatment planning station CRS-Master®, it offers an optimally coordinated system solution for refractive laser surgery – for a convenient workflow, efficient patient management and the best possible results.

The combination of precision and efficiency

- High-precision flaps due to high-performance femtosecond technology
- High reproducibility and consistency of flap thickness
- Easy repositioning of the flap
- Optimum workflow due to perfectly coordinated system components
- Smooth, finely structured surface of the stromal bed
- Prevents unnecessarily raised IOP and thus temporary losses of vision and trauma, due to the anatomically curved contact glass and the non-scleral suction.

Pivotal patient supporting system – Optimum workflows for greater comfort

The shared use of the pivotable patient supporting system by the VisuMax and the MEL 80 saves the patient from having to move from one treatment location to another. The patient experiences the surgery as an integrated process. Unnecessary waiting periods are minimized and treatment efficiency increased.

MEL 80 – All you need for optimum results

All the parameters of this high-precision excimer laser are geared to increasing efficiency, achieving optimum treatment results and rapid visual recovery. Key factors here are the extremely high ablation speed, customized treatment planning with the optional CRS-Master, the high-performance eyetracker system and eye registration.


CRS-Master – For a truly individual treatment

As a state-of-the-art treatment planning tool, the CRS-Master transfers data from the wavefront diagnosis and corneal topography into the MEL 80. This additional incorporation of patient-related data makes it possible to create a complete, individual eye profile for Femto-LASIK or Laser Blended Vision that is perfectly tailored to the patient.



Femto-LASIK


Laser-in-situ-Keratotomy




The VisuMax femtosecond laser creates the flap.



The patient moves to the MEL 80 excimer laser.



The flap is manually opened and folded back to expose the deeper corneal layer (stroma) beneath.



The MEL 80 excimer laser ablates the pre-calculated corneal tissue point by point.



The upper corneal layer is then repositioned following the refractive correction.

Keratoplasty

For high-precision tissue grafts

With the Keratoplasty option, the VisuMax® quickly converts into a state-of-the-art workstation for corneal grafts. The excellent cutting quality and laser control enable smooth lamellar and circular incision areas for high-precision results.



Broad spectrum for keratoplasty

The Keratoplasty option for the VisuMax comprises functions specially developed for corneal grafts and their optimization. High-precision and rapid incision progress, combined with a high level of reproducibility, mean that the VisuMax can be used for the three most important corneal graft procedures:

- Penetrating Keratoplasty (PKP)
- Deep Anterior Lamellar Keratoplasty (DALK)
- Descemet's Stripping Endothelial Keratoplasty (DSEK)

Perfect tissue grafts

The unique adapter attached to the headrest of the patient supporting system serves as an ideal work platform for preparing the corneal donor graft as well as the recipient's cornea.



The practical adapter provides a robust and sterile work surface for preparing the corneal graft.

Specially developed contact glass (type KP)

The curved contour of the contact glass prevents any unnecessary compression of the corneal tissue. It is also compatible with most artificial anterior chambers.



Contact glass (type KP) for preparing the donor cornea

Precise results

- High-precision cutting quality in anterior lamellar and endothelial Keratoplasty
- Penetrating Keratoplasty with perfect fit for donor and recipient cornea
- Precisely predictable incision pattern for greater reliability when preparing thin grafts
- Small spot distance for excellent cutting quality and easy separation of the tissue
- Flexible adjustment of the cutting parameters on the VisuMax

Maximum efficiency for optimum workflows

- 500 kHz laser pulse frequency for faster, more precise treatments and shorter cutting times (typically less than 60 seconds), even for very deep cuts
- High-quality surgical microscope for all treatment phases
- Practical adapter for preparing the donor cornea
- Special contact glass (type KP) compatible with most artificial anterior chambers




View through the surgical microscope of the VisuMax: Separation of the cut lamellae from the recipient cornea

The option "Endothelial Keratoplasty" is not intended for sale in the United States.

Incision for ICR

For flexible access

When implanting intracorneal ring (ICR) segments, surgeons benefit from the unique advantages of the VisuMax® femtosecond laser. The Incision for ICR option offers the possibility, for the first time, to create even inclined cutting geometries and tunnel segments between 90° and 270°. The repeatedly proven femtosecond laser technology ensures not only high-precision cutting quality but also previously unattainable degrees of freedom when defining the tunnel parameters.



Freely variable cutting parameters: for the first time, an incision can even be created parallel to the posterior surface of the cornea.

Degrees of freedom redefined

The VisuMax allows the corneal tunnels for the implantation of intracorneal ring (ICR) segments to be prepared quickly and easily, in the precise position and depth required. The wide range of adjustable parameters and the various combinations thereof give surgeons a unique degree of flexibility.

Tailor-made segments

For the first time, it is also possible to create tunnels with an arc angle of less than 360° using a femtosecond laser: With segments of between 90° and 270° partial tunnels can be designed individually, with precision and a high degree of flexibility. This means that individual intracorneal ring segments can also be implanted, without creating unnecessary incisions in the corneal tissue. Thus, tunnels can also be created for mixed ring segments with different geometries.



New treatment alternatives offer individual benefits

- The option to select between 0, 1 or 2 trapezoidal access incisions facilitates flexible tunnel access
- Width and inclination of the tunnel can be freely defined and precisely adjusted to the individual corneal shape and the applied ring geometry
- Seamlessly integrated into the user interface of the VisuMax, the ICR option provides maximum ease of use
- Rapid and intuitive entry of the necessary parameters
- The possibility to save user-defined cutting geometries as reusable templates increases workflow efficiency
- Reliability enhancement by graphic visualization of the parameter selection and automatic consistency check of the input parameters
- The unique contact glass is modelled on the natural shape of the corneal surface and thus contributes to patient-friendly treatment methods.
- Excellent control of tunnel preparation and ICR implant insertion as well as complete video documentation using the high-quality ZEISS surgical microscope

All common ICR products are supported.

The option Incision for ICR is not intended for sale in the United States.



Technical data

Installation and operating instructions

VisuMax femtosecond laser system

System components	Patient supporting system, including platform Integrated uninterruptible power supply (UPS) Surgical microscope with additional slit illumination Integrated with digital recording video camera
--------------------------	---

Laser parameters	Wavelength	1043 nm
	Pulse duration	220–580 fs
	Laser pulse rate	500 kHz

Installation and set-up conditions

Weight	855 kg (including patient supporting system, platform, UPS)
Recommended space requirement	3.50 x 4.20 m (3.60 x 6.00 m in combination with MEL 80 Excimer Laser)

Electrical connection	100–240 V, 50/60 Hz, max. 16 A Separately fused circuit
------------------------------	--

Operating conditions

Room temperature	18 to 25 °C
-------------------------	-------------

Atmospheric humidity	30 to 70 %
-----------------------------	------------

Accessories	Single-use contact glasses Treatment Pack (sizes S / M / L and type KP) Keratoplasty adapter for patient supporting system
--------------------	---



Your local contact:

Argentina

Carl Zeiss Argentina S.A.
Calle Nahuel Huapi 4015 / 25
C1430 BCO Buenos Aires
Argentina
Phone: +54 11 45 45 66 61
bruzzi@zeiss.com.ar

Australia

Carl Zeiss Pty Ltd
Tenancy Office 4, Level 1
40-52 Talavera Road
North Ryde NSW 2113
Australia
Phone: +61 2 9020 1333
med@zeiss.com

Austria

Carl Zeiss GmbH
Laxenburger Str. 2
1100 Vienna
Austria
Phone: +43 1 79 51 80
austria@zeiss.org

Belgium

Carl Zeiss NV-SA
Ikarooslaan 49
1930 Zaventem
Belgium
Phone: +32 2 719 39 11
info@zeiss.be

Brazil

Carl Zeiss do Brasil Ltda.
Av. Nações Unidas, 21711
CEP04795-100 São Paulo
Brazil
Phone: +55 11 5693 5521
medbrasil@zeiss.org

Canada

Carl Zeiss Canada Ltd.
45 Valleybrook Drive
Toronto, ON M3B 2S6
Canada
Phone: +1 800 387 8037
micro@zeiss.com

China

Carl Zeiss Shanghai Co. Ltd.
1/f., Ke Yuan Building
11 Ri Yin Nan Road
Waigaoqiao Free Trade Zone
2005 Yang Gao Bei Road
Shanghai 200131
China
Phone: +86 21 5048 17 17
sro@zeiss.com.cn

Czech Republic

Carl Zeiss spol. s r.o.
Radlická 14/3201
150 00 Prague 5
Czech Republic
Phone: +420 233 101 221
zeiss@zeiss.cz

France

Carl Zeiss Meditec France SAS
60, route de Sartrouville
78230 Le Pecq
France
Phone: +33 1 34 80 21 00
med@zeiss.fr

Germany

Carl Zeiss Meditec VG mbH
Carl-Zeiss-Strasse 22
73446 Oberkochen
Germany
Phone: +49 7364 20 6000
vertrieb@meditec.zeiss.com
Surgical Ophthalmology:
Phone: +49 800 470 50 30
iol.order@meditec.zeiss.com

Hong Kong

Carl Zeiss Far East Co. Ltd.
Units 11-12, 25/F
Tower 2, Ever Gain Plaza
No. 88 Container Port Road
Kwai Chung
Hong Kong
Phone: +852 2332 0402
czfe@zeiss.com.hk

India

Carl Zeiss India Pvt. Ltd.
22, Kensington Road
Ulsoor
Bangalore 560 008
India
Phone: +91 80 2557 88 88
info@zeiss.co.in

Italy

Carl Zeiss S.p.A.
Viale delle Industrie 20
20020 Arese (Milan)
Italy
Phone: +39 02 93773 1
post@zeiss.it

Japan

Carl Zeiss Meditec Japan Co. Ltd.
Shinjuku Ku
Tokyo 160-0003
22 Honchio-Cho
Japan
Ophthalmic instruments:
Phone: +81 3 33 55 0331
medsales@zeiss.co.jp
Surgical instruments:
Phone: +81 3 33 55 0341
cmskoho@zeiss.co.jp

Malaysia

Carl Zeiss Sdn Bhd.
Lot2, Jalan 243/51 A
46100 Petaling Jaya
Selangor Darul Ehsan
Malaysia
Phone: +60 3 7877 50 58
malaysia@zeiss.com.sg

Mexico

Carl Zeiss de México S.A. de C.V.
Avenida Miguel Angel de Quevedo
496
04010 Mexico City
Mexico
Phone: +52 55 59 99 0200
cz-mexico@zeiss.org

Netherlands

Carl Zeiss B.V.
Trapezium 300
Postbus 310
3364 DL Sliedrecht
Netherlands
Phone: +31 184 43 34 00
info@zeiss.nl

New Zealand

Carl Zeiss NZ Ltd
15B Paramount Drive
P.O. Box 121 - 1001
Henderson, Auckland 0650
Ph: +64 9 838 5626
med@zeiss.com

Poland

Carl Zeiss sp. z o.o.
ul. Lopuszanska 32
02-220 Warsaw
Poland
Phone: +48 22 858 2343
medycyna@zeiss.pl

Singapore

Carl Zeiss Ptd. Ltd.
50 Kaki Bukit Place
Singapore 415926
Singapore
Phone: +65 6741 9600
info@zeiss.com.sg

South Africa

Carl Zeiss (Pty.) Ltd.
363 Oak Avenue
Ferndale
Randburg 2194
South Africa
Phone: +27 11 886 9510
info@zeiss.co.za

South Korea

Carl Zeiss Co. Ltd.
Seoul 121-828
Mapo-gu
141-1, Sangsu-dong
2F, BR Elitel Bldg.
South Korea
Phone: +82 2 3140 2600
korea@zeiss.co.kr

Spain

Carl Zeiss Meditec Iberia S.A.
Ronda de Poniente, 15
Tres Cantos
28760 Madrid
Spain
Phone: +34 91 203 37 00
info@zeiss.es

Sweden

Carl Zeiss AB
Tegeludsvaegen 76
10254 Stockholm
Sweden
Phone: +46 84 59 25 00
info@zeiss.se

Switzerland

Carl Zeiss AG
Feldbachstrasse 81
8714 Feldbach
Switzerland
Phone: +41 55 254 7534
med@zeiss.ch

Thailand

Carl Zeiss Thailand
Floor 8, Thosapol Land Building 2
230 Ratchadapisek Road
Huaykwang, Bangkok 10310
Thailand
Phone: +66 2 2 74 06 43
thailand@zeiss.com.sg

United Kingdom

Carl Zeiss Ltd.
15-20 Woodfield Road
Welwyn Garden City
Hertfordshire, AL7 1JQ
United Kingdom
Phone: +44 1707 871200
info@zeiss.co.uk

United States of America

Carl Zeiss Meditec, Inc.
5160 Hacienda Drive
Dublin, CA 94568
USA
Phone: +1 925 557 4100
info@meditec.zeiss.com



Manufacturer:

Carl Zeiss Meditec AG

Goeschwitzer Strasse 51-52

07745 Jena

Germany

www.meditec.zeiss.com/VisuMax

Publication No: 000000-1980-341 Printed in Germany CZ-11/2012
The contents of the brochure may differ from the current status of approval of the product in your country. Please contact our regional representative for more information.
Subject to change in design and scope of delivery and as a result of ongoing technical development. Printed on elemental chlorine-free bleached paper.
© 2012 by Carl Zeiss Meditec AG. All copyrights reserved.